

# Nitric Oxide Synthase Inhibitor Screening Assay Kit

Catalog # KA6216

Size 1 Kit

## Applications

### Result data

Result data

□

DPI titrations: iNOS from mouse was incubated with various concentrations of DPI. Each concentration of inhibitor contained 20 v/v% DMSO (final 2 v/v% in 50 uL reaction).

## Specification

<b>Product Description</b>	Nitric Oxide Synthase Inhibitor Screening Assay Kit is used for screening nitric oxide synthase inhibitor.
<b>Sample Volume</b>	5 uL
<b>Detection Method</b>	Colorimetric
<b>Regulatory Status</b>	For research use only (RUO)
<b>Storage Instruction</b>	Store components of the kit at 4°C or -20°C as described in the protocol.
<b>Note</b>	<p>Result data</p> <p>Result data</p> <p>DPI titrations: iNOS from mouse was incubated with various concentrations of DPI. Each concentration of inhibitor contained 20 v/v% DMSO (final 2 v/v% in 50 uL reaction).</p>

## Applications

- Functional Study