Protein Expression Kit II (GST-Fusion Protein)

Catalog # KA1438 Size 1 Kit

Specification	
Product Description	This kit produces GST-fusion proteins at high purities and high yields. The wheat germ included in thi s kit is a specialized high yield wheat germ extract prepared by removing endogenous proteins as no nspecifically adsorb to GST-tag affinity purification resin.
Pack Size	Automatic synthesis 6 samples (6-mL translation)/run or 24 samples (1.2-mL translation)/run Manual synthesis 6-well plates, 24-well plates or 96-well plates.
Required Instrument	Automatic synthesis: Protemist® DT II Manual synthesis: Incubator (transcription at 37°C, translation at 15-26°C)
Regulation Status	For research use only (RUO)
Applied Instrument	Protemist® DT II
Storage Instruction	Store the component at -20°C or -80°C as shown on the protocol.

Applications

• Purification