

FFPE FISH PreTreatment Kit 2

Catalog # KA2376 Size 1 Kit

Specification	
Product Description	Pre-treatment is an essential step for effective FISH on Formalin-Fixed Paraffin-Embedded (FFPE) t issue sections. This kit contains key reagents used to pretreat FFPE tissue sections from dewaxing with xylene through proteolytic digestion to help permeabilization of cell membranes to facilitate pene tration of fluorescence labeled probes. After proteolytic digestion, proceed with subsequent steps with the FISH assay protocol.
Regulation Status	For research use only (RUO)
Supplied Product	Kit content: 1. 10x Paraffin Pretreatment Solution , 500 mL 2. Protease (Pepsin) , 0.5 mL x 10 vials
Storage Instruction	After arrival, please store the 10x Paraffin Pretreatment Solution at room temperature (RT), and st ore the Protease (Pepsin) at -20°C.

Applications

• Fluorescent In Situ Hybridization (Cell)

Protocol Download